

S. aureus meticilin resistente (SARM).

Staphylococcus aureus

Antibiograma

S -> Sensible; R -> Resistente

Staphylococcus aureus			CMI
Resistencia			
inducible a clindami+			Pos
Cefoxitin Screen	+		Pos
Clindamicina	R		0,25
Daptomicina	S		0,5
Eritromicina	R		≥ 8
Fosfomicina	S		32
Gentamicina	S		$\leq 0,5$
Levofloxacino	I		0,25
Linezolid	S		2
Ceftarolina	S		0,5
Mupirocina	S		≤ 1
Oxacilina	R		≥ 4
Penicilina G	R		$\geq 0,5$
Rifampicina	S		$\leq 0,03$
Teicoplanina	S		$\leq 0,5$
Tigeciclina	S		$\leq 0,12$
Tobramicina	S		≤ 1
Trimetoprima/			
Sulfametoxazol	S		≤ 10
Vancomicina	S		$\leq 0,5$

